

Agilent ScreenTape Consumables and Reagents for the Agilent TapeStation Systems



Agilent offers a variety of ScreenTape consumables and reagents for the Agilent TapeStation systems:

- RNA and High Sensitivity RNA ScreenTape

For the analysis of total RNA samples from eukaryotic or prokaryotic origin providing quality, quantity and sizing information. The RNA integrity number equivalent (RIN^e) delivers an instant and objective evaluation of eukaryotic and prokaryotic total RNA degradation.

- DNA and High Sensitivity DNA ScreenTape

For the separation and analysis of DNA fragments and libraries up to 5,000 base pairs. Ideal for sample quality control (QC) in next-generation sequencing (NGS) workflows.

- Genomic DNA ScreenTape

For the separation and analysis of DNA samples up to 60,000 base pairs. It provides an automated numerical assessment of genomic DNA quality, the DNA integrity number (DIN), and is therefore the ideal QC tool for NGS and array comparative genomic hybridization (aCGH) workflows.

- Cell-free DNA ScreenTape

For the separation and analysis of cell-free DNA (cfDNA) samples including detection of high molecular weight DNA contaminations. The software automatically calculates the %cfDNA of the sample.

TapeStation RNA ScreenTape and Reagents

The RNA ScreenTape provides a fully automated, efficient and reliable RNA analysis for RNA characterization and quality assessment. The RNA integrity number equivalent (RIN^e) provides an instant and objective evaluation of total RNA degradation.

RNA ScreenTape		
• H	He	
• H	- He	
·	- He	
+	The	
	- He	
	lite	
	- Ite	
-	- He	
-	like	
-	- Bre	
	- Ba	
-	- Be	
	- Ba	
	- Ba	
	il.	
	- Des	
	11.0	

RNA ScreenTape ordering information

Part number	Description	Quantitiy
5067-5576	RNA ScreenTape. For analysis of total RNA down to a sensitivity of 5 ng/ μ L. Includes 7 ScreenTape devices.	For 112 samples
5067-5577	RNA ScreenTape sample buffer For analysis of total RNA down to a sensitivity of 5 ng/µL. Includes 600 µL sample buffer. Order with 5067-5576.	For 112 samples
5067-5578	RNA ScreenTape ladder. For the analysis of total RNA down to a sensitivity of 5 ng/μL. Includes 10 μL ladder. Order with 5067-5576 and 5067-5577.	
5067-5579	High Sensitivity RNA ScreenTape. For the high sensitivity analysis of total RNA down to 100 pg/µL. Includes 7 ScreenTape devices.	For 112 samples
5067-5580	High Sensitivity RNA ScreenTape sample buffer. For the high sensitivity analysis of total RNA down to 100 pg/µL. Includes 250 µL sample buffer. Order with 5067-5579.	For 112 samples
5067-5581	High Sensitivity RNA ScreenTape ladder. For the high sensitivity analysis of total RNA down to 100 pg/μL. Includes 10 μL ladder. Order with 5067-5579 and 5067-5580.	

RNA ScreenTape specifications

Analytical specifications RNA ScreenTape		High Sensitivity RNA ScreenTape
Sensitivity ¹	5 ng/µL	100 pg/µL
Quantitative precision ³	10 % CV	15 % CV
Quantitative accuracy ³	±20 %	±30 %
Quantitative range ³	25 – 500 ng/µL	500 – 10,000 pg/μL
RIN ^e functional range ²	25 – 500 ng/μL	1,000 − 25,000 pg/µL
Maximum buffer concentration	200 mM Tris	10 mM Tris
in sample	20 mM EDTA	1 mM EDTA
Physical specifications		
Analyzia tima	16 samples < 20 min	16 samples < 35 min
	96 samples < 95 min	96 samples < 180 min
Samples per consumable	16	16
Sample volume required	1 µL	2 μL
Kit stability 6 months		6 months
Kit size 112 samples		112 samples

¹ Signal-to-noise >3 (single peak)

² RIN^e – RNA integrity number equivalent

³ Applicable to eukaryotic total RNA as sample

TapeStation DNA ScreenTape and Reagents

The DNA ScreenTape assays for the Agilent TapeStation systems are ideal for sample QC of input genomic DNA and downstream analysis within the Next Generation Sequencing workflow. Select the sizing range appropriate for your application.

D1000 ScreenTape		
•H)	(Br	
• H	- Br	
• H .	- Ha	
•H	- Br	
• H	- Br	
• H	- He	
• E .	- He	
• H	- Br	
• I	- Br	
• E •	- Br	
• = [R	
•==	- Br	
•==	- Re	
•==	R	
•==	R	
•==	R	

DNA ScreenTape ordering information

Part number	Description	Quantity	
D1000 ScreenTape assay – 35 bp to 1,000 bp			
5067-5582	D1000 ScreenTape. For the analysis of DNA from 35 to 1,000 bp. Includes 7 ScreenTape devices.	For 112 samples	
5067-5583	D1000 reagents. For the analysis of DNA from 35 to 1,000 bp. Includes ladder and sample buffer. Order with 5067-5582.	For 112 samples	
5067-5586	D1000 ladder. For the analysis of DNA from 35 to 1,000 bp. Includes 17 μ L ladder.		
5067-5602	D1000 sample buffer. For the analysis of DNA from 35 to 1,000 bp. Includes 400 μ L sample buffer.		
5067-5584	High Sensitivity D1000 ScreenTape. For the high sensitivity analysis of DNA from 35 to 1,000 bp. Includes 7 ScreenTape devices.	For 112 samples	
5067-5585	High Sensitivity D1000 reagents. For the high sensitivity analysis of DNA from 35 to 1,000 bp. Includes ladder and sample buffer. Order with 5067-5584.	For 112 samples	
5067-5587	High Sensitivity D1000 ladder. For the high sensitivity analysis of DNA from 35 to 1,000 bp. Includes 30 µL ladder.		
5067-5603	High Sensitivity D1000 sample buffer For the high sensitivity analysis of DNA from 35 to 1,000 bp. Includes 300 µL sample buffer.		
D5000 ScreenTape assay – 100 bp to 5,000 bp			
5067-5588	D5000 ScreenTape. For the analysis of DNA from 100 to 5,000 bp. Includes 7 ScreenTape devices.	For 105 samples	
5067-5589	D5000 reagents. For the analysis of DNA from 100 to 5,000 bp. Includes ladder and sample buffer. Order with 5067-5588.	For 105 samples	
5067-5590	D5000 ladder. For the analysis of DNA from 100 to 5,000 bp. Includes 10 μ L ladder.		
5067-5592	High Sensitivity D5000 ScreenTape. For the high sensitivity analysis of DNA from 100 to 5,000 bp. Includes 7 ScreenTape devices.	For 105 samples	
5067-5593	High Sensitivity D5000 reagents. For the high sensitivity analysis of DNA from 100 to 5,000 bp. Includes ladder and sample buffer. Order with 5067-5592.	For 105 samples	
5067-5594	High Sensitivity D5000 ladder. For the high sensitivity analysis of DNA from 100 to 5,000 bp. Includes 20 µL ladder.		
Cell-free DNA ScreenTape assay – 50 bp to 800 bp			
5067-5630	Cell-free DNA ScreenTape. For the analysis of cell-free DNA from 50 to 800 bp. Includes 7 ScreenTape devices.	For 112 samples	
5067-5631	Cell-free DNA reagents. For the analysis of cell-free DNA from 50 to 800 bp. Includes ladder and sample buffer. Order with 5067-5630.	For 112 samples	
5067-5632	Cell-free DNA ladder. For the analysis of cell-free DNA from 50 to 800 bp. Includes 50 µL ladder and sample buffer.		
5067-5633	Cell-free DNA sample buffer. For the analysis of cell-free DNA from 50 to 800 bp. Includes 300 µL sample buffer.		
Genomic DNA ScreenTape assay – 200 bp to > 60,000 bp			
5067-5365	Genomic DNA ScreenTape. For the analysis of genomic DNA from 200 to > 60,000 bp. Includes 7 ScreenTape devices.	For 105 samples	
5067-5366	Genomic DNA reagents. For the analysis of genomic DNA from 200 to > 60,000 bp. Includes ladder and sample buffer. Order with 5067-5365.	For 105 samples	

DNA ScreenTape specifications

Analytical specifications	D1000 ScreenTape	High Sensitivity D1000 ScreenTape	D5000 ScreenTape	High Sensitivity D5000 ScreenTape
Sizing range	35 – 1,000 bp	35 – 1,000 bp	100 – 5,000 bp	100 – 5,000 bp
Typical resolution	35 – 300 bp: 15% 300 – 1,000 bp: 10%	35 – 300 bp: 15% 300 – 1,000 bp: 10%	400 – 5,000 bp: 15%	400 – 5,000 bp: 15%
Sensitivity ¹	0.1 ng/µL	5 pg/µL	0.1 ng/µL	5 pg/µL
Sizing precision	5 % CV ²	5 % CV ²	5 % CV ⁴	10 % CV ⁴
Sizing accuracy	±10 % ^{2,3}	±10 % ^{2,3}	±10 % ⁴	±15 % ⁴
Quantitative precision	0.1 – 1 ng/μL: 15 % CV 1 – 50 ng/μL: 10 % CV	15 % CV	0.1 – 1 ng/μL: 15% CV 1 – 50 ng/μL: 10% CV	15 % CV
Quantitative accuracy	±20 % ²	±20 % ²	±20 %	±25 %
Quantitative range	0.1 – 50 ng/µL	10 – 1,000 pg/µL	0.1 – 50 ng/µL	10 – 1,000 pg/µL
Maximum buffer concentration in sample	20 mM KCl 60 mM phosphate buffer 60 mM guanidine-HCl 240 mM NaCl 60 mM NaOAc	7 mM KCl 20 mM phosphate buffer 20 mM guanidine-HCl 80 mM NaCl 20 mM NaOAc	250 mM KCl 250 mM Tris-HCl 125 mM NaCl 50 mM NaOAc 25 mM MgCl ₂ 25 mM BSA 25 mM guanidine-HCl	25 mM KCl 25 mM Tris-HCl 12.5 mM NaCl 5 mM NaOAc 2.5 mM MgCl ₂ 2.5 mM BSA 2.5 mM guanidine-HCl
Physical specifications				
Analysis time	16 samples < 20 min 96 samples < 90 min	16 samples < 20 min 96 samples < 105 min	15 samples < 25 min 96 samples < 135 min	15 samples < 20 min 96 samples < 120 min
Samples per consumable	16	16	15	15
Sample volume required	1 µL	2 µL	1 μL	2 µL
Kit stability	6 months	6 months	6 months	6 months
Kit size	112 samples	112 samples	105 samples	105 samples

¹ Signal-to-noise >3 (single peak)
² Measured using one ladder per ScreenTape device
³ Sizing accuracy for analysis with electronic ladder: ±20 %
⁴ Determined using ladder as sample

DNA ScreenTape specifications continued

Analytical specifications	Genomic DNA ScreenTape	Cell-free DNA ScreenTape
Sizing range	200 to > 60,000 bp	50 – 800 bp
Sensitivity ¹	0.5 ng/µL	20 pg/µL
Sizing precision ²	200 – 15,000 bp: 15 % CV	10 % CV
Sizing accuracy ²	200 – 15,000 bp: ±15 %	±15 % ³
Quantitative precision	15 % CV	10 % CV ²
Quantitative accuracy	±20 %	±20 % ²
Quantitative range	10 – 100 ng/µL	100 – 4000 pg/µL
DIN functional range ⁴	5–300 ng/µL	-
%cfDNA functional range	-	100 – 5000 pg/µl
Maximum buffer concentration	10 mM MgCl ₂	25 mM NaCl
in sample	50 mM NaCl	25 mM KCl
	10 mM NaOAc	3 mM EDTA
	10% ethanol	0.1 % NaN ₃
	10% 2-propanol,	5 mM phosphate buffer
	1 μg/μL glycogen	10% ethanol
		10% 2-propanol
Physical specifications		
Analysis time	15 samples: < 25 min	16 samples < 25 min
	96 samples: < 140 min	96 samples < 150 min
Samples per consumable	15	16
Sample volume required	1 µL	2 µL
Kit stability	6 months	6 months
Kit size	105 samples	112 samples

¹ Signal/noise ratio >3 (single peak)
² Determined using the ladder as sample
³ Sizing accuracy for analysis with electronic ladder: ±20 %
⁴ DIN - DNA Integrity Number



www.agilent.com/genomics/tapestation

For Research Use Only. Not for use in diagnostic procedures. PR7000-7791

This information is subject to change without notice.

© Agilent Technologies, Inc. 2020-2024 Published in the USA, November 1, 2024 5994-0264EN

